

Check List before starting an XAFS scan:

- ❑ Beam at the ring - on or off;
- ❑ HXMA Wiggler - on or off;
- ❑ Shutters - open or close;
- ❑ JJ X-ray Slits - open properly or close;
- ❑ Sample at proper detecting position, or not;
- ❑ Dead time at 32Ge detector (if fluorescence mode) - proper or not;
- ❑ Gain setup at SR570 - proper or not;
- ❑ Offset setup at SR570 - properly or not;
- ❑ Dark current - recorded and subtracted or not;
- ❑ Scaler status - at continuous mode or not;
- ❑ XAFS data collection configuration setup - properly or not;
- ❑ Name of the ExpData file - correct or not;
- ❑ Data storage path - correct or not;
- ❑ BLGraph configuration - correct or not;
- ❑ Detector configuration - correct or not.